Kit Gram-Nicolle

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Reference: 364320-0000

Kit Gram-Nicolle

STAINING KIT FOR THE DIFFERENTIATION OF GRAM-POSITIVE FROM GRAM-NEGATIVE BACTERIA

Principe:

Sudan black stain is insoluble in aqueous media, but soluble in lipidic and alcoholic media.

This property is used in anatomopathology for staining smears and tissue sections. In buffered phosphate medium and in presence of Alcohol, sudanophilic granules are stained and appear dark blue, whereas cytoplasm of the cells remain colourless.

Kit description:

Bottle 1 : Carbolic Gentian Violet	1 x 240 ml
Bottle 2 : Lugol, PVP-stabilized solution	1 x 240 ml
Bottle 3 : Fast differentiator	1 x 240 ml
Bottle 4 : Ziehl Carbolic Fuchsin 1/10	1 x 240 ml

Number of tests possible: 230 to 300 slides

Specific material required but not provided:

95°-Ethanol

Recommendations and/or notes for use:

Adding of Polyvinylpyrrolidone (PVP) to Lugol Solution helps to avoid iodine migration and then provides a satisfactory stability of the ready-to-use plastic bottle packaged solutions.

Specimen preparation:

Specimen must be treated in accordance with procedures available in the laboratory and promulgated by national authorities.

Use:

For Professional Use Only. **IVD**

The collection and processing of chemical and biological waste must be conducted by Organic phase Ethyl acetate based and registered companies.

Storage Temperature: 15-25°C

Opened bottle shelf life : report to expiration date on the label.

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

CLARK G., *Staining procedures*, Williams & Wilkins, 4th ed., 1981, p. 377-379.

Protocol N° 85

Kit Gram-Nicolle


STAINING KIT FOR THE DIFFERENTIATION OF GRAM-POSITIVE FROM GRAM-NEGATIVE BACTERIA

Staining procedure:

Processing time : 3 minutes


STEP 1 : Place the slide on a stand with smears on top. Cover the slide with **Carbolic Gentian Violet** (bottle 1) for 1 to 5 minutes. Get rid of the stain and Rinse with tap water to remove excess of Carbolic Gentian Violet.

STEP 2 : Rinse and cover with **Lugol, PVP-stabilized solution** (bottle 2) for 30 seconds to 1 minute. Thoroughly rinse with tap water.

STEP 3 and 4 : Discolour with the **Fast differentiator** (bottle 3). Quickly rinse with tap water. Cover the slide with **Ziehl Carbolic Fuchsin 1/10** (bottle 4) for 1 minute.

STEP 5 and 6 : Quickly rinse with tap water and let dry the smear. Microscopic examination with a x100-immersion-objective.

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